SPECTRO-SEPT

3124

**CLEANER • DISINFECTANT • DETERGENT
DEODORIZER • FUNGICIDE • MILDEWSTAT**

SPECTRO-SEPT is recommended for:

- hospitals;
- nursing homes;
- seniors residences;
- research and analysis laboratories;
- ambulance services;
- pharmaceutical industry;
- veterinary hospitals;
- kennels;
- school establishments;
- hotels;
- office buildings;
- establishments where housekeeping is of prime importance in controlling cross-contamination from treated surfaces.

ADVANTAGES

- Broad spectrum of activity.
- Outstanding cleaning efficiency.
- Neutral pH. Will not dull gloss floor finishes.
- Versatile. **SPECTRO-SEPT** is the only product needed for cleaning all surfaces that are washable with water.
- Rinse-free. Cleans and disinfects in one step.
- Effectively controls bad odours by killing microorganisms that cause bad odours while leaving a pleasant, fresh scent.
- Quaternary amines from 5th generation.
- Contains no phosphates, no solvents and no nonylphenols.
- Does not emit VOC.
- Biodegradable surfactants according to tests of the OECD 301 series.
- Effective bactericide and virucide in water up to 400 ppm hardness.

SPECTRO-SEPT is a phosphate-free neutral disinfectant cleaner designed to provide effective cleaning, deodorizing, and disinfection. This product, when used as directed, is formulated to disinfect hard non-porous inanimate environmental surfaces such as floors, walls, metal surfaces, bathrooms, shower stalls, bathtubs and cabinets. For larger areas such as operating rooms and patient care facilities, this product is designed to provide both general cleaning and disinfecting. This product deodorizes those areas which generally are hard to keep fresh smelling, such as garbage cans, storage areas, empty garbage bins and cans, toilet bowls, and other areas which are prone to odours caused by microorganisms.

DIRECTIONS

Application methods

Disinfection: To disinfect inanimate, hard non-porous surfaces add 8 ml of **SPECTRO-SEPT** per litre of water. Apply solution with a cloth, flip top dispensing cap, manual sprayer and wet mop so as to wet all surfaces thoroughly. Allow 10 minutes of contact time, then remove excess liquid. For heavily soiled areas, a pre-cleaning step is required. Always use a fresh solution.

To disinfect toilet bowls: Remove gross filth or soils from surfaces with bowl brush. Add 8 ml of **SPECTRO-SEPT** to the bowl water. Brush or swab the bowl completely using a scrub brush or toilet mop making sure to get under the rim. Let stand for 10 minutes and flush.

Deodorization: To deodorize, apply **SPECTRO-SEPT** as indicated under the heading DISINFECTION.

Bactericidal activity: At the 8 ml per litre dilution, **SPECTRO-SEPT** demonstrates effective activity against several organisms (see table on the following page).

Fungicidal activity: At the 8 ml per litre dilution, **SPECTRO-SEPT** demonstrates effective fungicidal activity against the pathogenic fungi Trichophyton mentagrophytes and Candida albicans.

Efficacy tests have demonstrated that this product is an effective fungicide in water up to 200 ppm hardness (as CaCO₃) in the presence of organic soil (5% blood serum).

Mildewstat: To control mold and mildew (such as Aspergillus niger) and the odours they cause on pre-cleaned, hard, non-porous surfaces add 8 ml per litre of water. Apply solution with a cloth, flip top dispensing cap, manual sprayer and wet mop making sure to wet all surfaces completely. Let air dry. Prepare a fresh solution for each use. Repeat application at weekly intervals or when mildew growth reappears.

SPECTRO-SEPT when used on environmental, inanimate hard non-porous surfaces at 8 ml per litre of water exhibits effective **virucidal activity** against several viruses (see table on the following page).

Kills HIV-1 and HIV-2 on pre-cleaned environmental surfaces/objects previously soiled with blood/body fluids in health care settings (Hospitals, Nursing Homes) or other settings in which there is an expected likelihood of soiling inanimate surfaces/objects with blood or body fluids, and in which the surfaces/objects likely to be soiled with blood or body fluids can be associated with the potential transmission of Human Immunodeficiency Virus Type 1 and Type 2 (HIV-1 and HIV-2).

Special instructions for cleaning and decontamination of HIV-1 and HIV-2 on surfaces soiled with blood or

Proportioning methods

body fluids.

Personal Protection: When handling items soiled with blood or body fluids, use disposable latex gloves, gowns, masks, or eye coverings.

Cleaning procedures: Blood or body fluids must be thoroughly cleaned from surfaces and objects before application of this product.

Contact Time/Dilution: At a use dilution of 8 ml of **SPECTRO-SEPT** per litre of water this product is effective against HIV-1 and HIV-2 in the presence of 5% blood serum when allowed a 10 minute contact time.

Disposal of Infectious Materials: Blood or other body fluids should be autoclaved and disposed of according to local regulations for infectious waste disposal.

Efficacy tests have demonstrated that this product is an effective bactericide and virucide in water up to 400 ppm hardness (as CaCO₃) in the presence of organic soil (5% blood serum).

PACKAGING

Case of 4 x 1 litre (33.8 US fl. oz.).
Case of 4 x 3.8 litres (1 US gallon).

BENEFITS OF THE DOSING BOTTLE

- Light, portable dosing system, easy to store on a cleaning cart.
- Graduated dosing from 10 to 40 ml.
- Easy and safe to use.
- Standardizes product usage by measuring equal doses time after time. Eliminates the use of over-concentrated solutions and reduces costs.
- Does not require an automatic dosing system and a dedicated water supply.

The information herein must not be reproduced without written permission. It represents our best judgment in the matter, but we can not assume responsibility for any loss or accident that may result from its use.

2013/08/23

**Products manufactured by: CHOISY LABORATORIES LTD., Louiseville (Québec) Canada J5V 2L7
CANADA: 819 228-5564 / 1 800 567-2157 • USA: 1 800 224-6061 • www.choisy.com • info@choisy.com**

CHEMISTRY AND BIOTECHNOLOGY SERVING SANITATION FOR A CLEAN ENVIRONMENT

SPECTRO-SEPT

3124

CHEMICAL AND PHYSICAL PROPERTIES

pH as is 25 °C: 6.0 - 7.0

Specific gravity 25 °C: 0.995 - 1.000

Viscosity 25 °C: less than 25 cP

Stability: stable at normal warehousing conditions

Physical properties: transparent yellow liquid

Odour: freshness

Active ingredients:

3.255% Octyl decyl dimethyl ammonium chloride

1.628% Diocetyl dimethyl ammonium chloride

1.628% Didecyl dimethyl ammonium chloride

4.339% Alkyl (40% C₁₂, 50% C₁₄, 10% C₁₆) dimethyl benzyl ammonium chloride

LEGISLATION

Exempt from WHMIS. The product holds a drug identification number (D.I.N.: 02295806).

Material safety data sheet available.

BIOCIDAL EFFICIENCY DATA

Bactericidal activity

At the 8 ml per litre dilution, **SPECTRO-SEPT** demonstrates effective activity against the following organisms:

ORGANISMS	ATCC#
Bordetella bronchiseptica	31437
Corynebacterium ammoniagenes	6871
Enterobacter aerogenes	13048
Enterobacter cloacae	23355
Enterobacter cloacae	clinical isolate
Enterococcus faecalis	19433
Enterococcus faecalis	clinical isolate
Enterococcus faecalis (VRE)	51575
Escherichia coli	11229
Escherichia coli	clinical isolate
Fusobacterium necrophorum	27852
Klebsiella pneumoniae subsp. pneumoniae	13883
Lactobacillus casei subsp. rhamnosus	7469
Listeria monocytogenes	35152
Pasteurella multocida	7707
Proteus mirabilis	25933
Proteus mirabilis	9921
Proteus vulgaris	13315
Pseudomonas aeruginosa	15442

ORGANISMS	ATCC#
Salmonella choleraesuis	10708
Salmonella choleraesuis subsp. choleraesuis serotype paratyphi B	8759
Salmonella choleraesuis subsp. choleraesuis serotype pullorum	9120
Salmonella choleraesuis subsp. choleraesuis serotype typhi	6539
Salmonella choleraesuis subsp. choleraesuis serotype typhimurium	14028
Serratia marcescens	8100
Shigella dysenteriae	12180
Shigella flexneri Type 2b	12022
Shigella sonnei	25931
Staphylococcus aureus	
Staphylococcus aureus	clinical isolate
Staphylococcus aureus (MRSA)	33592
Staphylococcus aureus (Vancomycin Intermediate)	
Staphylococcus aureus subsp. aureus	6538
Staphylococcus epidermidis	29461
Staphylococcus epidermidis	clinical isolate
Streptococcus pyogenes	clinical flesh eating strain - BIRD M3
Streptococcus pyogenes Group A	19615
Xanthamonas maltophilia	clinical isolate

Fungicidal activity

At the 8 ml per litre dilution, **SPECTRO-SEPT** demonstrates effective fungicidal activity against the following pathogenic fungi:

- Trycophyton-mentagrophytes
- Candida albicans

Mildewstat

At the 8 ml per litre dilution, **SPECTRO-SEPT** controls mold and mildew (such as Aspergillus niger).

Virucidal activity

At the 8 ml per litre dilution, **SPECTRO-SEPT** exhibits effective virucidal activity against several viruses:

- Bovine Rhinotracheitis
- Canine Distemper virus
- Feline Leukemia
- Feline Picornavirus
- Herpes Simplex Type 1 (causative agent for fever blisters)
- Herpes Simplex Type 2 (genital disease)
- Human Immunodeficiency Virus Type 1 and Type 2 (HIV-1 and HIV-2)
- Influenza A2/Hong Kong
- Pseudorabies
- Vaccinia

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